

## Pico Legend Specification

Model		PICO LEGEND
Wavelength		1064nm / 532nm
Operating Mode		585nm / 650nm / 595nm / 660nm
Pulse Energy	1064nm	Max. 500mj
	532nm	Max. 250mj
	585/595nm	Max. 150mj
	650/660nm	Max. 100mj
Pulse Width	Nd : YAG	350ps
Spot Size		2~10mm (1064 / 532nm) 2mm-3mm (585/650nm/595/660)
Pulse Rate		2Hz-10Hz 585 / 650 : 1-2 Hz
Aiming Beam		Diode 650nm
Optional Handpiece		· Dye Hand Piece (585,650nm / 595,660nm) · Fractional Handpiece



# Pico Legend

REAL  
PICO  
LASER

**01 Stable Operating**  
 ■ Minimal pain

**02 Convenient to Treat**  
 ■ Smart protocols

**03 Low Risk Of Side Effect**  
 ■ Fast treatment

**Before and After**

Scar Acne Scar		→	
Tattoo Removal		→	
Melasma		→	
Pigmen- tation		→	
Freckle		→	

■ **Advantage**

- Stable Operation
- Convenient to Treat with Smart Protocols
- Low risk of Side Effect
- Fast Treatment

■ **Pico Legend has the shortest pulse width of 350 ps**

■ **Low risk of side effect**

Less thermal effect on the treated area makes the best result of treatment with low risk of side effect


■ **Indications**

Tattoo, Wrinkle, Acne, Nevus of ota, Blue nevus, Speckled nevus, Lentigo, Cafe au lait, Hyperpigmentation, Freckle, Melasma and so on

■ **Super ultra short pulsed width**

Picosecond Nd:YAG Laser (1064nm&532nm + 595nm 660nm HP+ Fractional HP) with super ultra short pulse width

■ **HandPiece**



The fractional hand-piece uses fractional beams (360 spots, 16x16mm) to allow higher fluence (J/cm2) compared to the bulk mode hand-piece, allowing it to effectively destroy various pigmentations which were hard to treat before.



Fractional Hand Piece (1064nm)    Zoom Hand Piece (2-10mm Standard)    Dye Hand Piece (660 / 650 / 595 / 585nm)

